

SpeckNet Ltd	IF1.33 10 Crichton Str EH8 9AB Edinburgh	tel 01316505182
		jmann@inf.ed.ac.uk

Respeck System Design

Prepared for: Edwin Lindsay, Principal Consultant

Prepared by: Janek Mann, Research Associate

26 January 2012


SpeckNet Ltd	IF1.33 10 Crichton Str EH8 9AB Edinburgh	tel 01316505182
		jmann@inf.ed.ac.uk

System Design

Objective

The Respeck wireless respiratory monitoring system is design to monitor respiratory rate and activity in a non-obtrusive fasion.



Overall System



The Respeck system consists of one or a small number of Respeck wireless patches connected via 2.4GHz radio to a basestation. The basestation can either operate independently powered from a standard USB charger (microUSB connector) sending data to remote servers via GPRS or storing to microSD card, or with a USB connection to a Home Health Hub or PC.

SpeckNet Ltd	IF1.33 10 Crichton Str EH8 9AB Edinburgh	tel 01316505182
		jmann@inf.ed.ac.uk


Respeck Patch



The Respeck Patch consists of an ARM microcontroller, a non-rechargeable Lithium polymer battery, an accelerometer and a low-power 2.4GHz radio.

SpeckNet Ltd	IF1.33 10 Crichton Str EH8 9AB Edinburgh	tel 01316505182
		jmann@inf.ed.ac.uk

Basestation



The Respeck basestation consists of an ARM core microcontroller, power management unit with battery charger, microUSB socket for charging and data exchange, microSD card interface (non user serviceable), SIM card socket (non user serviceable), GSM/GPRS modem, 2.4GHz low-power radio, RF PA and LNA (combined package), button and indicator lights.

Two antennas are used (for GSM bands and 2.4GHz). Connection to the antennas is via u.FL connectors.